



# Reclining Wheelchair

## User Manual



**FOR MODELS# 11950, 11952**

**BEFORE** using this product, read this manual and save for future reference.

## **Table of Contents**

<b>Wheelchair Diagram.....</b>	<b>1</b>
<b>Safety Precaution.....</b>	<b>2</b>
<b>Opening Wheelchair.....</b>	<b>3</b>
<b>Folding Wheelchair.....</b>	<b>4</b>
<b>Removing Arms.....</b>	<b>4</b>
<b>Swing Back A4 Swing Back Arms.....</b>	<b>5</b>
<b>Removing Front Riggings.....</b>	<b>5</b>
<b>Adjusting Footplates.....</b>	<b>6</b>
<b>Adjusting Leg Rests.....</b>	<b>6</b>
<b>Applying Wheel Locks.....</b>	<b>7</b>
<b>Reclining Back Section.....</b>	<b>8</b>

# Wheelchair Diagram



## **Safety Precautions**

Safety is of paramount importance in all activities outlined in this manual. It is crucial to thoroughly familiarize yourself with secure methods when participating in routine tasks. Attendants must also be well-versed in and adhere to proper assistance protocols.

For personalized approaches tailored to your specific requirements, we strongly recommend consulting your physician, nurse, or physical therapist.

Please note that wheelchairs are not designed for use as seats in motor vehicles.

**Engaging Wheel Locks** - Prior to sitting down or rising from the chair, ensure that the wheel locks are engaged against both larger wheels.

**Ensuring Stability** - Under no circumstances should you shift your position on the seat by leaning forward outside of the chair. Should it become necessary to do so, ensure that the casters are positioned forward. This action elongates the wheelbase, enhancing the wheelchair's stability. Additionally, refrain from reaching beyond the length of your arm, as this action may result in tipping.

**Transferring** - Avoid stepping on the footplates while standing up, as this action can lead to tipping. Either fold up the footplates or move them to the side.

**Navigating Ramps** - Prior to attempting an ascent on a ramp, it's important to understand your personal physical capabilities and limitations in terms of strength and endurance. To master proper techniques for ramp negotiation, it's advisable to seek guidance from a certified professional.

**Dealing with Steep Inclines:** In situations involving inclines exceeding 10% (equivalent to a 1-foot elevation on a 10-foot ramp), it's essential to have an attendant's assistance to prevent loss of steering or tipping backward.

During incline ascents, if you need to halt, exercise caution and avoid sudden or abrupt forward movements when resuming motion, as such actions can lead to the wheelchair tipping over.

When going down an incline, grasp the handrims as you descend to maintain speed and control. DO NOT use the wheel locks to slow down for this may cause accidental locking, sudden stopping, or even tipping to one side. Also, try and avoid changing directions while descending as this can cause instability.

## **Opening Wheelchair**

Grab the armrest with your one hand and tilt the wheelchair towards you. Push down the seat rail with your other hand. \*Do not put your fingers between the seat rail and the side panel.



## Folding Wheelchair

If footrests or leg rests are on the wheelchair, fold the foot- plates before folding the wheelchair. Grab the armrest with your one hand and tilt the wheelchair towards you. Lift up the front of the seat rail opposite you with your other hand.



## Removing Arms

To remove arm, press down the button lock with your finger & lift the armrest.



## Swing Back Arms

To swing back arms, press the lever to release the lock and swing back the arms.



## Removing Front Riggings

To remove the footrest or leg rest, push the lock handle towards the inside of the wheelchair. To lock the footrest or leg rest, swing the footrest or leg rest towards the inside of the wheelchair.



## Adjusting Footplates

Fold up the footplate and loosen the hex nut; after adjusting the height, tighten the hex nut.



## Adjusting Leg Rests

To elevate leg rests:

Simply tilt the leg rest with your hand to the desired position.





To lower leg rests:

Push the release lever to lower the leg rest. hold the leg rest with your other hand to prevent sudden dropping.



## **Applying Wheel-Locks**

Push the handle forward until it locks. Pull the handle back to release the wheel locks.



## **Reclining Back Section**

Squeeze both handles at the same time and recline the backrest to desired angle. Release handles when the backrest reaches the desired position.



## ***LIMITED WARRANTY***

Your product is warranted to be free of defects in materials and workmanship for three years of the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment, This 3 year Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service,

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips. which are subject to normal wear and need periodic replacement.