



Serenity

Alternating Pressure Low Air Loss Mattress System Cell on Cell

The Selectis Serenity alternating low air loss mattress system helps provide prevention therapy and treatment of multiple stage wounds. Cell on Cell design prevents “bottoming out” while ensuring mattress remains inflated in case of power failure.



SPECS

Dimensions	61758	35" x 76" x 8"
	61058	35" x 80" x 8"
	61458	35" x 84" x 8"
Weight Capacity	350 lbs	
Warranty	18 Months	

FEATURES

- Cell on Cell design ensures that mattress remains inflated
- 8" deep air bladders provide supreme comfort & pressure redistribution
- Adjustable setting allows for customized pressure by resident weight
- 10-minute cycle for alternating pressure
- 18 laser air holes provide a cool and dry environment
- Audio alarm signals low pressure in the mattress
- Static mode for easy resident transfer
- Quilted 2 way stretch fluid resistant, vapor permeable cover provides extra comfort
- CPR valve ensures rapid deflation
- Quick release tubing from the pump