

Patient Lift Slings

Residents who need assistance with mobility will benefit from the lift sling, which is compatible with other manufacturer's lifts. With a variety of styles and models, caregivers can select the type that best suits their resident's needs.

Mesh Full Body Sling

Full body slings are ideal for residents who need full support, especially for the head and neck. The mesh material offers more breathability and dries easily when exposed to moisture.

Mesh Full Body Sling with Commode

Full body sling with commode option provides easy toileting access, and can also be used for residents who need to be in a constant seated position.

Solid Full Body Sling

Polyester full body slings are more durable and provide full body support. This sling is ideal for transfer from a wheelchair to a bed or another chair, or transfer from floor to bed.



Patient Lift Slings

Solid Divided Leg Sling

The divided leg sling is the most versatile option, and is quick and easy to use. It helps the residents distribute their weight more comfortably.



Sit to Stand Sling

The Sit to Stand sling is used for transfer from a seated position in and out of bed to a chair, wheelchair, toilet, or from chair to chair. It is designed for residents with minimal weight-bearing ability, who have head and neck control, can bend at the hips, knees and ankles, and can sit up on the edge of the bed.



Sit to Stand Buttock Strap

The Sit-to-Stand Buttock Strap is a supportive accessory designed to complement our sit-to-stand sling. It enhances stability during transfers and offers adjustable comfort. Pair it with our sling for safer and easier mobility assistance.



Spec Table

* See next page for reference

| Model # | Type | Size | Binding Color | Dimension A (in) * | Dimension B (in) * | Dimension C (in) * | Weight Capacity |
|---------|------------------------------------|------|---------------|--------------------|--------------------|--------------------|-----------------|
| #71236 | Full Body Sling, Mesh | M | Purple | 39 | 51 | 23 | 600 lbs |
| #71237 | | L | Green | 43 | 55 | 26 | 600 lbs |
| #71238 | | XL | Blue | 42 | 58 | 26 | 600 lbs |
| #71239 | | XXL | Black | 44 | 52 | 44 | 600 lbs |
| #71256 | Full Body Sling with Commode, Mesh | M | Purple | 39 | 51 | 24 | 600 lbs |
| #71257 | | L | Green | 43 | 55 | 24 | 600 lbs |
| #71258 | | XL | Blue | 43 | 59 | 27 | 600 lbs |
| #71259 | | XXL | Black | 43 | 51 | 44 | 600 lbs |
| #71241 | Full Body Sling, Solid | M | Purple | 40 | 51 | 26 | 600 lbs |
| #71242 | | L | Green | 43 | 56 | 27 | 600 lbs |
| #71243 | | XL | Blue | 43 | 60 | 26 | 600 lbs |
| #71244 | | XXL | Black | 46 | 64 | 28 | 600 lbs |
| #72101 | Divided Leg Sling, Solid | S | Navy | 25 | 57 | 35 | 600 lbs |
| #72102 | | M | Purple | 22 | 59 | 39 | 600 lbs |
| #72103 | | L | Green | 26 | 62 | 43 | 600 lbs |
| #72104 | | XL | Blue | 29 | 69 | 46 | 600 lbs |
| #72105 | | XXL | Black | 32 | 72 | 49 | 600 lbs |
| #72109 | Sit to Stand Sling | M | N/A | 16 | 32 | N/A | 600 lbs |
| #72110 | | L | N/A | 16 | 44 | N/A | 600 lbs |
| #72112 | | XL | N/A | 16 | 57 | N/A | 600 lbs |
| #72113 | | XXL | N/A | 19 | 63 | N/A | 600 lbs |
| #72115 | Buttock Strap | N/A | N/A | 7.5 | 34.5 | N/A | 400 lbs |




Sit to Stand
Sling



Full Body
Sling, Solid



Full Body
Sling, Mesh



Divided Leg
Sling



Full Body
Sling with
Commode,
Mesh



Buttock
Strap

